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Intensiv rotierende diamantbeschichtete Dentalinstrumente

Produktbeschreibung

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Table with columns for ISO size (1/10 mm), speed range, and contact force (Max, N) for various grit sizes (Körnung Extra fine, Fine, Medium, Coarse).

Instructions for use

Intensiv rotating diamond-coated instruments for Dentistry

Product description

- Stainless steel instruments with diamond coated working part
Diamond grit-sizes: 8µm – 150µm
Shapes: ball, inverted cone, pear, cylinder, flame, tapered, pointed, torpedo, football, wheel, lenticular, special shapes
Available in FG (Friction Grip) and RA (Right Angle)
Compatible with turbines and/or contra-angles
Sterilizable and reusable

Indications

Abrasion of dental material for reducing, contouring, finishing and/or polishing in Dentistry

Application details

- Preparations and Finishing for restorations and prosthetics in Dentistry
Root surface treatments in Periodontics

Instructions for use

- Insert the instrument as deeply as possible into the chuck (FG) or until safety "click" is noticed (RA) (figure 1a).
Set the handpiece at the required working speed (table 1) before applying the instrument on the treated area.
Activate water spray during the entire treatment (minimum 50 mL/min) for cooling and instruments cleaning.
Ensure that there is consistent water flow, regulate in case of water excess or shortage.
Provide additional water spray cooling for instruments exceeding working part diameter ISO 020 and working part length of 10 mm, and diamond grain-size smaller or equal to 50 microns.
Use the instruments applying a contact force up to the maximum values listed in table 1.
It is recommended to use a dental dam during treatment.
Whenever clinically possible, use the entire length of the working part to allow for lower contact pressure and a uniform usage of the instrument (figure 1b).
Once the preparation has been completed, remove the instrument from the treated area and bring it to a standstill.
It is recommended to wear gloves and safety glasses.

Maintenance and sterilization

- Pre-clean by disinfection the instruments immediately after each usage.
Clean the instruments and remove possible debris after each use (with soft brush, cleaning rubber Intensiv Diakleen or ultrasonic bath, to maintain their abrasive properties.
Instruments are packed non-sterile. They must be disinfected and sterilized prior to first use on the patient and disinfected, cleaned with soft brush or sonic bath and sterilized immediately after each use.
Disinfect the diamond-coated instruments separately from other instruments such as polymeric polishers and abrasives stones.
Use cleaning and disinfection solutions strictly observing the concentrations and reaction times indicated by the manufacturer.
In case of the presence of grinding residues on the instruments it is advisable to use an ultrasonic bath, cleaning rubber Intensiv Diakleen and/or soft brush for cleaning.
After disinfection and cleaning, inspect the instruments for residues. If necessary, repeat the disinfection/cleaning procedure.

- Check for possible damage; dispose of oxidized, eccentric or deformed instruments, and instruments where the diamond-coated working part is totally or partially worn-out.
Sterilization must be carried out according to validated procedures.
Cleaning and disinfection can be also carried out using professional thermosdisinfectors.
Concerning the sterilization process, refer to ISO standard 17664.

Risk warnings

- Avoid jamming or levering actions when rotating, as this increases the risk of instrument breakage.
Never exceed the specified maximum speed, to avoid instrument breakage caused by the generation of powerful centrifugal forces.
Partial consumption of the diamond-coated working part can cause localized overheating during usage of these instrument.
Avoid applying forces greater than the recommended values, as this could cause damage to the instrument and the treated area.
The use of instruments with fine grit size could lead to overheating of the tooth.
To avoid loss of traceability of the instruments, during their entire application it is necessary to keep the packaging with indicated lot number.

Figure 1


Table 1. Speed and contact force for Intensiv rotating diamond-coated instruments

Table with columns for ISO size (1/10 mm), speed range, and contact force (Max, N) for various grit sizes (Grit Extra fine, Fine, Medium, Coarse).

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