## Oscillating Instruments

DE

### Oszillierende Instrumente

Intensiv Bevelshape  
Intensiv Cavishape  
Intensiv Proxoshape  
Intensiv Proxoshape Coarse  
Intensiv Proxoshape Flexible  
Intensiv Rootshape

Zu verwenden in Kombination mit:

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To be used in combination with:

- Intensiv Swingle, WG-69 LT Cons, with light
- Intensiv Swingle, WG-69 A Cons, without light

Alternatively:

- KaVo 61 LRG or 61 LG
- KaVo contra-angle E61C

#### Intended use

Abrasion of dental material for reducing, contouring, finishing and/or polishing in Dentistry

#### Application details

- Intensiv Bevelshape:  
Margin beveling in cavity, crown and chamfer preparations  
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C
- Intensiv Cavishape:  
Contouring and finishing of cavity margins  
Intensiv Swingle WG-69 LT / KaVo E61C / KaVo 61LG
- Intensiv Proxoshape:  
Contouring, finishing and pre-polishing of proximal filling surfaces  
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C
- Intensiv Proxoshape Coarse:  
Modelling of interproximal aesthetic restorations  
Intensiv Swingle WG-69 LT / KaVo E61C / KaVo 61LG
- Intensiv Proxoshape Flexible:  
Finishing and polishing of aesthetic restorations with difficult to access proximal areas. Disposal of excess and over hangings filling  
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C
- Intensiv Rootshape:  
Rootplaning and odontoplasty  
Intensiv Swingle WG-69 LT / KaVo 61 LRG / KaVo 61 LG / KaVo E61C

#### Recommendations for use

- Application by dentist only.
- Set the instrument to the recommended speed before applying it to the tooth.
- Do not stop the instrument in the interdental space.
- In case of very tight contact points: separate the interdental space by placing a wooden wedge cervically.
- Insert the instrument as far as possible.
- Set the micro motor at 40'000 RPM max. With Intensiv Swingle this corresponds to 20'000 oscillations per minute approx. (2:1).
- Provide an adequate water spray cooling (minimum 50ml/min.).
- Regulate water flow in case of water excess or shortage.
- Pay attention to applied load (around 50 grams), up to a maximum bending of the strip of approximately 1 mm.
- Remove the strip from the interdental space whilst the contraangle is moving.

#### Maintenance and Sterilization

- Clean the instruments and remove debris after each use, so as to maintain their abrasive properties.
- Instruments are delivered non sterile. They must be sterilized prior to first use and disinfected, cleaned, and sterilized after each use.
- Use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Hot air sterilizers are not suitable for diamond instruments.
- Sterilization must be carried out according to validated procedures.
- Full cleaning, disinfection, sterilization, checking and maintenance instructions can be requested at info@intensiv.ch

#### Risk Warnings

- Avoid bending and levering over the limit: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Oscillating over 20'000 oscillations per minute could cause breakage of the instrument.
- Avoid applying excessive loads that would lead to the bending of the instrument, as this could cause damage to the instrument itself.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging.

In the unlikely case of an incident caused by the device please refer to the manufacturer and competent authorities.

Intensiv SA declines any liability for any use of the instrument other than that specified in these instructions for use.

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