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Instructions for use

Mode d'emploi

DE

Intensiv ProxoshapeTray

Product description

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Bur-shank silicone holders orange, Ref. 056orange/10
Bur-shank silicone holders yellow, Ref. 056yellow/10
Bur-shank silicone holders white, Ref. 056white/10
Bur-shank silicone holders red, Ref. 056red/10
Bur-shank silicone holders brown, Ref. 056brown/10
Bur-shank silicone holders grey, Ref. 056grey/10
Bur-shank silicone holders dark blue, Ref. 056darkblue/10
Bur-shank silicone holders green, Ref. 056green/10
Bur-shank silicone holders light blue, Ref. 056lightblue/10
Bur-shank silicone holders assorted, Ref. 056ass/8

GB

Intensiv ProxoshapeTray

Product description

- Tray in stainless steel, consisting in a base and lid connected by a hinge.
- The Tray can be closed by overlapping the cover on the base.
- Once it is closed, the instruments inside are stable in their position and will not move.
- The Tray has an internal folding support with nine holes.
- The holes are suitable to support sterilisable coloured silicone holders.
- The silicone holders enable to contain the instruments in a compact and steady arrangement, suitable for the safe and orderly withdrawal and relocation of the instruments.
- The colour of the silicone holders can be chosen to match the colour code of the instrument diamond grit-size.
- Dimensions:
 - Length 87mm
 - Depth 51mm
 - Height 14mm
- The Tray has slots for the disinfection liquid to flow through it.
- The Tray is designed to hold, through its silicone holders, nine Intensiv oscillating files, such as: Intensiv Proxoshape, Intensiv Proxoshape Coarse, Intensiv Proxoshape Flexible, Intensiv Bevelshape, Intensiv Rootshape, Intensiv Cavishape

Indications

- Arrangement of Intensiv oscillating files during interproximal preparation, excess removal, finishing and polishing procedures.
- Storage of Intensiv oscillating files.

Practical tips

- To open the Tray lift the cover upwards, slightly pressing the top wing and holding firmly the base.
- Move the internal folding support upwards, in order to have it parallel to the working surface. Two little wings on left and right side of the support help to tilt it in that direction.
- Place the silicone holders inside the nine holes before inserting the instruments, as shown in the figures.
- Make sure that the folding internal support is moved back downwards before the Tray is closed.
- Close the Tray by pressing the cover downwards until you hear a "click".

Maintenance and Sterilisation

- Intensiv ProxoshapeTray is delivered non-sterile. It must be disinfected and sterilized prior to first application to patient.
- It must be disinfected, cleaned and sterilized immediately after each use.
- If disinfection by immersion is carried out, it is appropriate to only use cleaning/disinfection solutions that provide protection against corrosion, carefully following the amounts indicated by the manufacturer.
- Sterilisation must take place according to established procedures. If possible, use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Chemiclave sterilisers may also be used.
- For sterilisation procedures refer to the ISO 17664 standard; the following suggestions are recommended:
 - Cycles at 134°C:
 - Tmin = 134°C – Tmax= 138°C
 - Pressure = 3.15 bar abs
 - Time = 4 min (can be incremented)
 - Cycles at 121°C:
 - Tmin = 121°C – Tmax = 125°C
 - Pressure = 2.10 bar abs
 - Time = 16 min (can be incremented)

Intensiv ProxoshapeTray, Ref. PST500

A wide variety of different colour bur-shank silicone holders is available:

Bur-shank silicone holders orange, Ref. 056orange/10
Bur-shank silicone holders yellow, Ref. 056yellow/10
Bur-shank silicone holders white, Ref. 056white/10
Bur-shank silicone holders red, Ref. 056red/10
Bur-shank silicone holders brown, Ref. 056brown/10
Bur-shank silicone holders grey, Ref. 056grey/10
Bur-shank silicone holders dark blue, Ref. 056darkblue/10
Bur-shank silicone holders green, Ref. 056green/10
Bur-shank silicone holders light blue, Ref. 056lightblue/10
Bur-shank silicone holders assorted, Ref. 056ass/8

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Bur-shank silicone holders orange, Ref. 056orange/10
 Bur-shank silicone holders yellow, Ref. 056yellow/10
 Bur-shank silicone holders white, Ref. 056white/10
 Bur-shank silicone holders red, Ref. 056red/10
 Bur-shank silicone holders brown, Ref. 056brown/10
 Bur-shank silicone holders grey, Ref. 056grey/10
 Bur-shank silicone holders dark blue, Ref. 056darkblue/10
 Bur-shank silicone holders green, Ref. 056green/10
 Bur-shank silicone holders light blue, Ref. 056lightblue/10
 Bur-shank silicone holders assorted, Ref. 056ass/8

Bur-shank silicone holders positioning and removal