DE

## Intensiv ProxoStrip Anterior Pat. 699 819

### Produktbeschreibung

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## Intensiv ProxoStrip Anterior Pat. 699 819

### Product description

- Metal strip, diamond-coated on one side only with two working parts of different grit size: 40µm (red) and 15µm (yellow)
- Perforated anatomical ends for a better gripping
- Height: 2.5mm
- Total length: 55mm
- Thickness: 0.05mm
- Sterilizable and reusable

### Indications

- Finishing and polishing of interproximal aesthetic restorations for anterior teeth
- Beveling of the interproximal areas

### Practical Tips

- The reduced length of the strip allows an optimal finger support on the anterior adjacent teeth.
- To get past the contact area, use the uncoated part of the strip and exert a slight horizontal movement.
- After introduction, work with the less abrasive side (yellow, 15µm), in order to facilitate the access of the more abrasive side (40µm, red). Then begin the finishing procedure using a back-and-forward movement.
- Next, switch again to the less abrasive side (15µm, yellow) to complete the procedure with the polishing.
- ProxoStrip Anterior must be removed with a vertical movement after repositioning it to the uncoated segment.
- In the event the contact area is very tight, it is useful to start by inserting a wooden wedge to open up a slight separation between the teeth before placing ProxoStrip Anterior in the area to be treated.

### Maintenance and sterilisation

- Intensiv ProxoStrip Anterior is supplied non-sterile. Prior to first use on the patient and immediately after each use, the instruments must be disinfected, cleaned, and sterilized.
- Keep Intensiv ProxoStrip Anterior separate from non-stainless steel instruments, such as rubber polishers and abrasive stones.
- If disinfection by immersion is being used, it is important to use only cleaning/disinfection solutions that provide corrosion protection and strictly observe the concentrations recommended by the manufacturer.
- Sterilization must be performed with validated procedures. If possible, use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Chemiclave sterilizers may also be used.
- Avoid temperatures above 180°C which may affect the durability of the instruments.
- Concerning the sterilization process we refer to the ISO standard 17664; so we suggest to follow these advices:
  - Cycles at 134°C:  
Tmin = 134°C – Tmax = 138°C  
Pressure = 3.15 bar abs  
Time = 4 min (raisable)
  - Cycles at 121°C:  
Tmin = 121°C – Tmax = 125°C  
Pressure = 2.10 bar abs  
Time = 16 min (raisable)

### Intensiv ProxoStrip Anterior (grit size 40µm and 15µm)

Ref. PXA4015/6 Package of 6 Strips  
Ref. PXA4015/12 Package of 12 Strips

### Additional models

#### Intensiv ProxoStrip (grit size 40µm and 15µm)

Ref. PX4015/6 Package of 6 Strips  
Ref. PX4015/12 Package of 12 Strips

#### Intensiv ProxoStrip Plus (grit size 40µm and 15µm)

Ref. PXP4015/6 Package of 6 Strips  
Ref. PXP4015/12 Package of 12 Strips

#### Intensiv ProxoContour (grit size 60µm and 40µm)

Ref. PXC6040/6 Package of 6 Strips

#### Intensiv ProxoContour Coarse (grit size 80µm and 40µm)

Ref. PXC8040/6 Package of 6 Strips

#### Intensiv ProxoPolish (grit size 15µm and 8µm)

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